

Improve Stem Cell Yield

Sci-tive-N

Stem Cell Investigations - Total In Vitro Environment

Sci-tive has been designed to mimic in-vivo conditions providing a continuous cell culture environment which eliminates cellular stress linked to variations in temperature, pH and oxidation yielding better stem cell lines.

🌀 “Lab in a Box” Concept

Reduces need for laminar air flow cabinets and CO₂ incubators

🌀 Stress-Free Cell Manipulation

Now Stem Cell Biologists can work under stable O₂, CO₂, pH, temperature and humidity for all manipulations, mixing and culturing steps

🌀 Direct-Hand Access into Workstation

Using sterilized and user friendly cuff & sleeve system (Ezee Sleeve™)

🌀 Aids Compliance with EU TCD (2004/23/EU)

VOC filtration and HEPA filtration create a biologically safe working environment

🌀 CE (1293) Approved (MDD 93/42/EEC)

Compliant with European EMC and electrical safety tests



innovation in incubation

Product Specifications

External Workstation Dimensions

- Height of Unit (excluding stand) 40"
- Width of Unit (side to side) 65"
- Depth 30"
(additional 6" required behind workstation for piping)

Internal Workstation Dimensions

- Height 23"
- Width 46"
- Depth 22"

Workstation Weight

- Weight approx. 500 lbs

Interlock Dimensions

- Height 8"
- Width 10"
- Depth 7"

Atmosphere Control

- Temperature Ambient + 5°C - 45°C
- Humidity Control Ambient - 85%
- O₂ Control 0.0% - 20.9%
- CO₂ Control 0.0% - 30.0%
- Integrated HEPA Filtration System
- Activated Carbon Filtration System

Standard Accessories

Workstation

- Stand
- 2 Racks with Perforated Shelves to Optimise Homogeneity
- Direct-hand Access using Ezee Sleeve™ and Ezee Pleez System
- Removable Front with UV Filtration Screen
- Waste Trays
- Gas Tight Waste Port

Gas Control Software

- USB Communications Port
- 7 Days Continuous Storage of Event Log Data
- One Touch Calibration (O₂)
- On-Screen Fault Assistance


Alarms

- Temperature - Visual and Audible
- Humidity - Visual
- Gas Low Pressure - Visual and Audible

Sci-tive-N

Stem Cell Investigations - Total In Vitro Environment

The Ruskinn brand was founded in 1993 and rapidly became established as one of the world's leading suppliers and manufacturers of anaerobic and modified atmosphere workstations. Ruskinn's purpose built factory is located within the Sony Technology Centre near Pencoed in South Wales, UK. Ruskinn's range of high quality workstations is available worldwide through a global network of carefully selected distributors.



NEW Optional Accessories

- Remote Monitoring

BAKER BIO SCIENCE SOLUTIONS

Sold exclusively in the United States by Baker BioScience Solutions
A Division of The Baker Company, Inc.

Baker BioScience Solutions
The Baker Company, Inc.
P.O. Drawer E
161 Gatehouse Road
Sanford, Maine 04073 USA
sales@bakerbioscience.com
www.bakerbioscience.com

 **RUSKINN**
TECHNOLOGY LTD